



Setup video

# Serene Air

True low air loss alternating pressure therapy & pulsation system

Serene Air is a dual mode system with Alternating Pressure Therapy (APT) or Pulsation options. Serene Air also has a Constant Low Pressure (static) option, as well as a max inflate and seat inflation for the Fowler position. It includes a zippered full-length pocket with a 1" high foam base. Initial inflation takes less than three (3) minutes, and the surface has a quilted top cover with high moisture vapor transmission rate (MVTR). Serene Air has several added safety and convenience features: auto-lock on the settings after five (5) minutes, accessible and highly visible CPR release valve at the control pump, plus a handy carrying bag for compact storage and transport.

## Features at a glance

- Suitable for vulnerable clients at high to very high risk for pressure injuries
- Provides effective microclimate control through low air loss technology
- High performance blower offers >120 L/min of air loss for effective microclimate control
- Fluid-resistant, quilted cover is made from Silver+ nylon, enhancing durability and hygiene
- Zippered cover allows for one-step removal
- Support clients weighing up to 920 lb (specifically on the 48-inch mattress)
- Available in three widths: 36 inches, 42 inches, and 48 inches, with a length of 79 inches and a height of 8 inches
- Available sizes:
  - 36" / 42" / 48" x 79" L x 8" H
- Optional gap filler extends the surface to a length of 84"



Control panel