



Prevention Through
Innovation.

PRODUCT FEATURES

- 3" wide transfer edges enables safe and easy transfers on and off the surface reducing falls and entrapment.
- Suitable for clients with high to very high risk for pressure injuries.
- Available in 36" / 42" / 48" widths.

SERENE **EDGE** True Low Air Loss System

The Serene *EDGE* dual mode system with Alternating Pressure Therapy or Pulsation features a high performance blower that offers 120 lpm of Low Air Loss for an effective microclimate control.

ADVANCED PUMP TECHNOLOGY

Multiple modes include Alternating Pressure Therapy and Constant Low Pressure (static). Serene *EDGE* also features max inflate and seat inflation for the Fowler position. Client weight references are incorporated into the control panel so caregivers can easily find the correct level of comfort and firmness needed for optimum comfort and healing.

ACTIVE THERAPY IMPROVES HEALING

Serene *EDGE* Air Pulsation and Alternating Therapy Pressure (APT) modes assist in wound healing and prevention. The APT has four operating cycles (5/10/15/20min), while pulsation increases and decreases the surface's air volume every 15 seconds.

SAFETY AND CONVENIENCE

To protect your clients, the Serene *EDGE* has a number of added safety features. An accessible and highly visible CPR release valve allows caregivers to quickly deflate the surface during emergencies. Plus there is an auto-lock on the settings after five (5) minutes. When setting it up, initial inflation takes less than three (3) minutes. The surface also has a quilted top cover with high moisture vapor transmission rate.



Size: 80" Length x 36" Wide x 7" High
Weight Limit: 450 lbs.

Warranty: Mattress: 36 months, Cover: 6 months

**Contact your medical equipment supplier
for more information or to purchase.**

DISTRIBUTED BY:

5 MINDS MOBILITY GROUP

ggarland@5mindsmobility.com

orders@5mindsmobility.com

555 Washington Ave., Suite 310

St. Louis, MO 63101

Ph: 866-709-1332 Fax: (877)-426-5697